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PZN 02486923 4209002

Glucosefine® Glucose Test Strips

IMPORTANT: Please read this information and your Glucosefine® operation manual before using the Glucosefine® glucose test strips.

Intended Use:

The Glucosefine® Glucose Test Strips are used with Glucosefine® Glucose Test Meter for the measurement of glucose in your blood. The Glucosefine® Blood Glucose Monitoring System is plasma-calibrated to allow comparison of results with the laboratory method. The Glucosefine® Glucose Test Strips are for testing outside the body (in vitro diagnostic use only). Glucosefine® Blood Glucose monitoring system is for self-testing.

General :

Glucosefine® glucose test strips adopt a new bio-sensor technology composed of electro-chemical components and require a small 0.5µl volume of blood giving you less pain. Glucosefine® can store 365 readings in memory with date & time, you can easily store your readings on your PC (optional). This will provide an easier and better way to control your glucose levels.

Storage and Handling :

- Store the Glucosefine® Glucose Test Strip vials in a cool and dry place between 2-30°C(36-86°F). Keep out of direct sunlight. Do not freeze.
- Store Glucose Test Strips in their original vial only. Do not mix the Glucose Test Strips in new vials or in any other container.
- Immediately replace the cap and close tightly after removing a Glucosefine® Glucose Test Strip.
- Use each Glucose Test Strip immediately after removing it from the vial.
- Do not use Glucose Test Strips after the expiry date printed on the package or vial since it may cause inaccurate results.
- Make a notation of the discard date on the vial label when you first open it.
- Discard remaining Glucosefine® Glucose Test Strips 3 months after first opening the vial.
- Avoid getting dirt, food, and water on the Glucose Test Strip. Do not handle Glucose Test Strips with wet hands.
- Do not bend, cut, or alter any Glucosefine® Glucose Test Strip.
- Glucosefine® Glucose Test Strips are for single use only. Do not re-use.
- Do not perform tests at temperatures below 10°C(50°F) or above 40°C(104 °F). Do not perform tests with humidity below 10 % or above 90 %.
- Be careful not to pollute the color coding label. (Polluted color coding label due to dust and other substances)
- Discard the used Glucose Test Strip carefully, to prevent any infection.
- Keep the Test Strips away from children!

| |
|---|
| Warning! To avoid false readings, use only a Glucosefine® Meter to measure glucose with your Glucosefine® Glucose Test Strips. |
|---|

Care Procedures :

- Do not reuse glucose test strips. Single use only.
- If you experience any symptoms that are not consistent with your blood glucose test results and you have followed the instructions described in your Glucosefine® Meter Manual, call your physician.
- Do not make significant changes to your diabetes control program without consulting your physician.
- Do not ignore physical symptoms without consulting your physician.

Test Principle :

Glucose in the blood sample will react to the electrodes in the glucose test strip, generating an electrical current that will stimulate a chemical reaction. This reaction is measured by the Glucosefine® meter and displayed as your blood glucose result.

Reagent Composition :

Each Glucosefine® glucose test strip contains:

| | |
|--|-----------|
| Glucose Oxidase (Aspergillus niger) : | 1.5 Units |
| Mediator (Hexaaminmeruthenium chloride) : | 100 µg |
| Binder : | 3.16 µg |
| Stabilizer : | 47.7 µg |

The procedure for Blood Glucose Measurement :

Materials needed: Glucosefine® Test Strips, Glucosefine® meter, operation manual, lancing device, sterile lancet

Obtaining Blood Sample :

The Glucosefine® Glucose Test Strips are designed to be used with fresh capillary whole blood. You may obtain a blood drop from either a fingertip or an alternative site (see meter manual).

To obtain a drop of blood, follow these steps :

- Step 1: Cleanliness** Wash your hands and test site with warm water. Dry it thoroughly. (Dirt, oil, lotion, or sweat may affect the test result).
- Step 2: Lancing** Prepare the lancing device. Insert a lancet into the lancing device. The device holds, positions, and controls how deeply the lancet goes into the skin. (Refer to the owner’s manual for more information).
- Step 3: Insert the Glucose Test Strip** Remove a test strip from the vial. Recap the vial immediately to prevent moisture from affecting the other strips. Insert the test strip into the port of the meter, with arrow sign facing up and showing towards the meter. The meter will automatically turn on and display the code number. Make sure that the code on your meter display matches the code number printed on the test strip vial. If not matching repeat this step. The test strip icon with a drop of blood will be flashing at the top of the meter display. You are now ready for testing.
- Step 4: Obtain and apply the blood sample** Place the lancing device against the area to be pricked, for example the finger tip. For a correct test result, an amount of 0,5µl is necessary. While the symbol for apply is flashing, place the drop of blood on your finger against the slim end of the test strip until the window is completely filled. Do not press your finger against the test strip. Do not try to wipe off the blood with the test strip. A beep sound confirms, when the window is correctly filled. If you have a poor amount of blood: Run a new test with a NEW test strip!

Note: Pricking the side of your fingertip is less painful. Do not squeeze or milk the puncture site.

Step 5: Measurement result in 5 seconds

The meter automatically starts testing. It counts down from 5 to 1, displays the measurement result and stores it automatically (To store with User’s Activity and User-ID please refer to your meter manual). Use the test strip ejector to remove the test strip. The meter will automatically turn off.

Important Information About Using Alternate Site Testing :

- Alternate sites where you can test are dorsal hand, ventral palm, upper arm, forearm, calf, and thigh.
- Under certain conditions, blood glucose test results obtained using samples taken from your alternate sites may differ significantly from fingertip samples.
- When blood glucose is changing rapidly, fingertip samples show these changes more quickly than alternate sites samples.
- When your blood glucose is falling, testing with a fingertip sample may identify a hypoglycemic (low blood sugar) level sooner than a test with a alternate sites sample.
- Use alternate sites samples only for testing prior to or more than two hours after meal, insulin doses, or physical exercise.
- Testing performed within two hours after meals, insulin doses, or physical exercise or whenever you feel that your glucose levels may be changing rapidly should be done from the fingertip.
- You should also use fingertip testing whenever you have a concern about hypoglycemia (insulin reactions) such as when driving a car, particularly if you suffer from hypoglycemic unawareness (lack of symptom to indicate as insulin reaction), as forearm testing may fail to detect hypoglycemia.

Test Results :

The results are displayed on the Glucosefine® meter as either milligrams of glucose per deciliter of blood (mg/dL) or millimoles of glucose per liter of blood (mmol/L). The glucose test meter displays results between 10–600mg/dL (0,6–33,3mmol/L).

If the test result is below 10mg/dL(0,6mmol/L), "Lo" will appear on the Glucosefine® meter display. Consult with your doctor immediately in how to treat extremely low blood glucose (hypoglycemia). If the test result is above 600mg/dL(33,3mmol/L), "Hi" will appear on the Glucosefine® meter display. This indicates extremely high blood glucose (hyperglycemia). Seek medical assistance immediately.

IMPORTANT NOTE: Blood glucose may be altered by your food, physical activity, and/or stress.

Range of Expected Values :

Self testing of blood glucose levels provides a way to control your diabetes. Consult with your physician to determine the best range of expected blood glucose values for you.

The fasting adult blood glucose range for a person without diabetes is less than 100mg/dL or 5.5 mmol/L. Two hours after meals, normal glucose levels should be less than 140mg/dL or 7.8mmol/L.

* American Diabetes Association: Diabetes Care, Volume 34, Supplement 1, January 2011, p. S11- S61.

IMPORTANT:

If you have a test result below 60mg/dL(3,3mmol/L) or over 240mg/dL(13,3mmol/L) contact your doctor immediately. If your blood glucose result is unusually low or high, repeat the test again with a new glucose test strip. If the results are still inconsistent, please consult your physician before making any decision to control your diabetes.

Quality Control (System Maintenance) :

Glucosefine® Glucose Control Solution is used to check that the meter and glucose test strips are properly working together to provide you with the most accurate reading.

Note:

The Glucosefine® Glucose Control Solutions is sold separately. The normal, and high level glucose control solutions can be obtained through your pharmacy, Metrado GmbH or your local representative.

Glucosefine® Glucose Control Solution test can only be used with the Glucosefine® meter and should be used during the following:

- When a new vial of glucose test strips is opened.
- Any suspicion that the meter or glucose test strips are not working properly.
- When your blood glucose test results are not consistent with your symptoms, or if you think they are not accurate.
- If you drop the glucose test meter.
- Use routinely to obtain accurate result.

When the glucose control solution is applied to the top edge of the glucose test strip, you should get results within the expected range printed on the label of the glucose test vial. If the glucose control solution test results fall outside the range, repeat the

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